

Anti-Aif1/Iba1 Polyclonal Antibody For Wb, If-P, Elisa - P55008

Item Number: CM0016597



Introduction

Rabbit polyclonal antibody against human AIF1/IBA1, validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for macrophage activation and vascular inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	AIF1/IBA1 Rabbit pAb
Remarks/Alias	IBA1; IRT1; AIF-1; IRT-1; AIF1/IBA1
Species	Human
Gene ID (Human)	199
Gene ID	199
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human AIF1/IBA1 (NP_001614.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55008
Protein Weight	17kDa
Shipping	Ice bag