

Anti-Aif1/Iba1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P55008

Item Number: CM0020763



Introduction

Human AIF1/IBA1 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Ideal for macrophage activation, vascular inflammation, and lymphocyte migration research studies.

[Learn More](#)

Parameter	Value
Product Name	AIF1/IBA1 Rabbit pAb
Remarks/Alias	IBA1; IRT1; AIF-1; IRT-1; AIF1/IBA1
Species	Human
GeneID (Human)	199
GeneID	199
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human AIF1/IBA1 (NP_001614.3))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P55008
Protein Weight	17kDa
Shipping	Ice bag