

Anti-Aifm2/Fsp1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Brq8

Item Number: CM0006512



Introduction

Anti-AIFM2/FSP1 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets ferroptosis suppressor protein 1 (FSP1), a key inhibitor of ferroptosis. Molecular weight 41kDa.

[Learn More](#)

Parameter	Value
Product Name	AIFM2/FSP1 Rabbit mAb
Remarks/Alias	AMID; FSP1; PRG3; AIFM2/FSP1
Species	Human
GeneID (Human)	84883
GeneID	84883
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-371 of human AIFM2/FSP1 (NP_116186.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRQ8
Protein Weight	41kDa
Shipping	Ice bag