

Anti-Aifm2/Fsp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Brq8

Item Number: CM0016848



Introduction

High-quality rabbit polyclonal antibody targeting human AIFM2/FSP1 (ferroptosis suppressor protein 1), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for ferroptosis and redox biology research.

[Learn More](#)

Parameter	Value
Product Name	AIFM2/FSP1 Rabbit pAb
Remarks/Alias	AMID; FSP1; PRG3; AIFM2/FSP1
Species	Human
GeneID	84883
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-371 of human AIFM2/FSP1/AMID/AMID (NP_116186.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRQ8
Protein Weight	41kDa
Shipping	Ice bag