

Anti-Aim2 Rabbit Polyclonal Antibody For Wb, Elisa - Q91Vj1

Item Number: CM0023590



Introduction

Anti-Aim2 rabbit polyclonal antibody validated for Western blot and ELISA. Detects mouse Aim2 and cross-reacts with rat. Suitable for studying inflammasome activation, DNA damage response, and tumor suppression mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Aim2 Rabbit pAb
Remarks/Alias	Gm1313; Ifi210; Aim2
Species	Mouse
Gene ID	383619
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of mouse Aim2 (NP_001013801.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q91VJ1
Protein Weight	40kDa
Shipping	Ice bag