


## Anti-Aim2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O14862

Item Number: CM0025816



### Introduction

High-quality rabbit polyclonal antibody against human AIM2, validated for WB, IF/ICC, and ELISA. Detects AIM2 in human, mouse, and rat samples. Ideal for inflammasome and innate immunity research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | AIM2 Rabbit pAb                                                                                                           |
| Remarks/Alias    | PYHIN4; AIM2                                                                                                              |
| Species          | Human                                                                                                                     |
| GeneID (human)   | 9447                                                                                                                      |
| GeneID           | 9447                                                                                                                      |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 130-337 of human AIM2 (NP_004824.1). |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, IF/ICC, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | O14862                                                                                                                    |
| Protein Weight   | 39kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |