


Anti-Aimp1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12904

Item Number: CM0027083



Introduction

Rabbit polyclonal antibody targeting human AIMP1, suitable for WB, IHC-P, and ELISA. Detects human, mouse, and rat AIMP1. Ideal for studies on angiogenesis, inflammation, and tRNA synthetase biology.

[Learn More](#)

Parameter	Value
Product Name	AIMP1 Rabbit pAb
Remarks/Alias	p43; HLD3; EMAP2; SCYE1; EMAPII; AIMP1
Species	Human
Gene ID (Human)	9255
Gene ID	9255
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-180 of human AIMP1 (NP_001135888.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12904
Protein Weight	34kDa
Shipping	Ice bag