

Anti-Aip Rabbit Polyclonal Antibody For Wb And Elisa - 000170

Item Number: CM0017277



Introduction

Rabbit polyclonal antibody targeting human AIP (ARA9/XAP2), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying AHR signaling and HBV X protein regulation.

[Learn More](#)

Parameter	Value
Product Name	AIP Rabbit pAb
Remarks/Alias	ARA9, XAP2, PITA1, XAP-2, FKBP16, FKBP37, SMTPHN, AIP
Species	Human
Gene ID (Human)	9049
Gene ID	9049
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human AIP (NP_003968.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00170
Protein Weight	38kDa
Shipping	Ice bag