

Anti-Aipl1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nzn9

Item Number: CM0029660



Introduction

Rabbit polyclonal antibody targeting human AIPL1 (aryl-hydrocarbon-interacting protein-like 1), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for vision research and Leber congenital amaurosis studies.

[Learn More](#)

Parameter	Value
Product Name	AIPL1 Rabbit pAb
Remarks/Alias	LCA4; AIPL2; AIPL1
Species	Human
GeneID (Human)	23746
GeneID	23746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-384 of human AIPL1 (NP_055151.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NZN9
Protein Weight	44kDa
Shipping	Ice bag