

Anti-Ajuba Rabbit Monoclonal Antibody (Kd Validated) For Wb, Ip, Elisa - Q96If1

Item Number: CM0006429

Introduction

High-quality rabbit monoclonal antibody against human Ajuba (JUB), KD validated for WB, IP, and ELISA. Detects human and rat Ajuba (57 kDa). Ideal for studying cell adhesion, cytoskeletal dynamics, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Ajuba Rabbit mAb
Remarks/Alias	JUB; [KD Validated] Ajuba
Species	Human
GeneID	84962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 201-290 of human Ajuba (NP_116265.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96IF1
Protein Weight	57kDa
Shipping	Ice bag