


## Anti-Ak1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00568

Item Number: CM0013632



### Introduction

Anti-AK1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA in human and mouse. AK1 is a cytoplasmic adenylate kinase involved in energy homeostasis and nucleotide metabolism.

[Learn More](#)

| Parameter        | Value                                                                                                                               |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | AK1 Rabbit pAb                                                                                                                      |
| Remarks/Alias    | HTL-S-58j; AK1                                                                                                                      |
| Species          | Human                                                                                                                               |
| GeneID (human)   | 203                                                                                                                                 |
| GeneID           | 203                                                                                                                                 |
| Immunogen        | Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human AK1 (NP_000467.1) |
| Source           | Rabbit                                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                                               |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse                                                                                                                        |
| SWISS            | P00568                                                                                                                              |
| Protein Weight   | 22kDa                                                                                                                               |
| Shipping         | Ice bag                                                                                                                             |