


## Anti-Ak1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P00568

Item Number: CM0014544



### Introduction

High-quality rabbit polyclonal antibody against human adenylate kinase 1 (AK1). Suitable for WB, IHC-P, ELISA in human, mouse, and rat samples. Ideal for cellular energy homeostasis and adenine nucleotide metabolism research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | AK1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | HTL-S-58j; AK1                                                                                                  |
| Species          | Human                                                                                                           |
| GeneID (human)   | 203                                                                                                             |
| GeneID           | 203                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human AK1 (NP_000467.1). |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IHC-P, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | P00568                                                                                                          |
| Protein Weight   | 22kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |