

Anti-Ak2 Polyclonal Antibody For Wb And Elisa - P54819

Item Number: CM0026304



Introduction

Rabbit polyclonal antibody against human AK2, validated for WB and ELISA. Useful for researching mitochondrial energy homeostasis, nucleotide metabolism, and hematopoiesis. Cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	AK2 Rabbit pAb
Remarks/Alias	ADK2; AK2
Species	Human
GeneID (Human)	204
GeneID	204
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-232 of human AK2 (NP_001616.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54819
Protein Weight	26kDa
Shipping	Ice bag