


## Anti-Ak3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uij7

Item Number: CM0029827



### Introduction

High-quality rabbit polyclonal antibody against human AK3, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for mitochondrial adenylate kinase research and metabolic studies.

[Learn More](#)

| Parameter        | Value                                                                                                          |
|------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name     | AK3 Rabbit pAb                                                                                                 |
| Remarks/Alias    | AK6; FIX; AK3L1; AKL3L; AKL3L1; AK3                                                                            |
| Species          | Human                                                                                                          |
| GeneID (Human)   | 50808                                                                                                          |
| GeneID           | 50808                                                                                                          |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human AK3 (NP_057366.2) |
| Source           | Rabbit                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                          |
| Application      | WB, ELISA                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                              |
| SWISS            | Q9UIJ7                                                                                                         |
| Protein Weight   | 26kDa                                                                                                          |
| Shipping         | Ice bag                                                                                                        |