


Anti-Ak4 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P27144

Item Number: CM0018329



Introduction

AK4 Rabbit Polyclonal Antibody validated for WB, IP, ELISA. Detects human and mouse mitochondrial adenylate kinase 4 (AK4), a key enzyme in cellular energy homeostasis. Ideal for metabolism research.

[Learn More](#)

Parameter	Value
Product Name	AK4 Rabbit pAb
Remarks/Alias	AK3; AK 4; AK3L1; AK3L2; AK4
Species	Human
GeneID (Human)	205
GeneID	205
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human AK4 (NP_037542.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P27144
Protein Weight	25kDa
Shipping	Ice bag