

## Anti-Akap14 (Akap28) Rabbit Monoclonal Antibody For Wb, Elisa - Q86Un6

Item Number: CM0003729

### Introduction

High-quality rabbit monoclonal antibody against human AKAP14 (AKAP28, PRKA14). Validated for WB and ELISA, suitable for studies on PKA anchoring and ciliary function. Reacts with human, mouse, rat.

[Learn More](#)

| Parameter        | Value                                                                          |
|------------------|--------------------------------------------------------------------------------|
| Product Name     | AKAP14 Rabbit mAb                                                              |
| Remarks/Alias    | AKAP28, PRKA14, AKAP14                                                         |
| Species          | Human                                                                          |
| GeneID (Human)   | 158798                                                                         |
| GeneID           | 158798                                                                         |
| Immunogen        | Synthetic peptide corresponding to amino acids 70-150 of human AKAP14 (Q86UN6) |
| Source           | Rabbit                                                                         |
| Category         | Monoclonal Antibodies                                                          |
| Application      | WB, ELISA                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                              |
| SWISS            | Q86UN6                                                                         |
| Protein Weight   | 23kDa                                                                          |
| Shipping         | Ice bag                                                                        |