

Anti-Akap4 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q5Jqc9

Item Number: CM0014955

Introduction

Rabbit polyclonal antibody against human AKAP4 (A-kinase anchor protein 4). Validated for WB, IF-P, ELISA; cross-reacts with human, mouse, rat. Detects a ~94 kDa testis-specific sperm fibrous sheath protein.

[Learn More](#)

Parameter	Value
Product Name	AKAP4 Rabbit pAb
Remarks/Alias	HI; p82; CT99; FSC1; PRKA4; AKAP-4; AKAP82; AKAP 82; hAKAP82; AKAP4
Species	Human
GeneID (Human)	8852
GeneID	8852
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 189-340 of human AKAP4 (NP_003877.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5JQC9
Protein Weight	94kDa
Shipping	Ice bag