

## Anti-Akap5 Rabbit Polyclonal Antibody For Wb, Elisa - P24588

Item Number: CM0014055



### Introduction

High-quality rabbit polyclonal anti-AKAP5 antibody validated for Western blot and ELISA. Recognizes human and mouse AKAP5, a scaffold protein in cAMP signaling and synaptic plasticity. Ideal for neuroscience research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | AKAP5 Rabbit pAb                                                                                                 |
| Remarks/Alias    | H21; AKAP75; AKAP79; AKAP5                                                                                       |
| Species          | Human                                                                                                            |
| GeneID (human)   | 9495                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human AKAP5 (NP_004848.3) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse                                                                                                     |
| SWISS            | P24588                                                                                                           |
| Protein Weight   | 47kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |