

Anti-Akap7 Rabbit Polyclonal Antibody For Wb, Elisa - O43687 / Q9P0M2

Item Number: CM0014882

Introduction

High-quality rabbit polyclonal antibody targeting human AKAP7 (A-kinase anchor protein 7), suitable for Western blotting and ELISA applications. Detects mouse and rat orthologs, aiding studies of PKA-mediated signaling. - O43687 / Q9P0M2

[Learn More](#)

Parameter	Value
Product Name	AKAP7 Rabbit pAb
Remarks/Alias	AKAP15, AKAP18, AKAP7
Species	Human
GeneID (human)	9465
GeneID	9465
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-81 of human AKAP7 (NP_004833.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O43687/Q9P0M2
Protein Weight	8kDa/11kDa/39kDa
Shipping	Ice bag