

Anti-Akap8 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43823

Item Number: CM0021837

Introduction

Rabbit polyclonal antibody against human AKAP8 (AKAP95). Suitable for WB, IHC-P, and ELISA; cross-reacts with human and rat. Ideal for chromatin condensation and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	AKAP8 Rabbit pAb
Remarks/Alias	AKAP-8; AKAP95; AKAP 95; AKAP-95; AKAP8
Species	Human
GeneID (Human)	10270
GeneID	10270
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 592-692 of human AKAP8 (NP_005849.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O43823
Protein Weight	76kDa
Shipping	Ice bag