

## Anti-Akr1A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P14550

Item Number: CM0014400



### Introduction

Human AKR1A1 polyclonal antibody validated for WB, IHC-P, ELISA. Detects aldo-keto reductase with roles in detoxification and carbonyl reduction. Reacts with human, mouse, rat.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | AKR1A1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | ALR; ARM; DD3; ALDR1; HEL-S-6; AKR1A1                                                                             |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 10327                                                                                                             |
| GeneID           | 10327                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human AKR1A1 (NP_006057.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | P14550                                                                                                            |
| Protein Weight   | 37kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |