

Anti-Akr1A1 Rabbit Polyclonal Antibody [Ko Validated] For Wb And Elisa - P14550

Item Number: CM0024241



Introduction

Rabbit polyclonal antibody targeting human AKR1A1, validated for WB and ELISA, cross-reactive with mouse and rat. Suitable for studying aldo-keto reductase function, detoxification, and S-nitrosylation regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] AKR1A1 Rabbit pAb
Remarks/Alias	ALR; ARM; DD3; ALDR1; HEL-S-6; A1
Species	Human
Gene ID	10327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human AKR1A1 (NP_006057.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14550
Protein Weight	37kDa
Shipping	Ice bag