

Anti-Akr1B1 Monoclonal Antibody For Wb, If/Icc, Elisa - P15121

Item Number: CM0006693



Introduction

Rabbit monoclonal antibody against AKR1B1 (aldose reductase) for WB, IF/ICC, ELISA. Detects human, mouse, rat AKR1B1 protein, ~36 kDa, suitable for metabolism and diabetic complication studies.

[Learn More](#)

Parameter	Value
Product Name	AKR1B1 Rabbit mAb
Remarks/Alias	AR; ADR; ALR2; ALDR1; AKR1B1
Species	Human
GeneID (Human)	231
GeneID	231
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human AKR1B1 (NP_001619.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15121
Protein Weight	36kDa
Shipping	Ice bag