

Anti-Akr1B1 Polyclonal Antibody For Wb And Elisa - P15121

Item Number: CM0015909



Introduction

Rabbit polyclonal antibody against human AKR1B1 (Aldose reductase). Validated for WB and ELISA, with cross-reactivity to mouse. Suitable for metabolism and diabetic complication studies.

[Learn More](#)

Parameter	Value
Product Name	AKR1B1 Rabbit pAb
Remarks/Alias	AR; ADR; ALR2; ALDR1; AKR1B1
Species	Human
GeneID (Human)	231
GeneID	231
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-316 of human AKR1B1 (NP_001619.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P15121
Protein Weight	36kDa
Shipping	Ice bag