


Anti-Akr1B10 Polyclonal Antibody For Wb, Ip, Elisa - O60218

Item Number: CM0027070



Introduction

Rabbit polyclonal antibody targeting AKR1B10 (ARL-1) for human, mouse, rat. Validated in WB, IP, ELISA. Detects aldo-keto reductase family 1 member B10 (36 kDa, O60218).

[Learn More](#)

Parameter	Value
Product Name	AKR1B10 Rabbit pAb
Remarks/Alias	HIS; HSI; ARL1; ARL-1; ALDRLn; AKR1B11; AKR1B12; AKR1B10
Species	Human
Gene ID (Human)	57016
Gene ID	57016
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-316 of human AKR1B10 (NP_064695.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60218
Protein Weight	36kDa
Shipping	Ice bag