


Anti-Akr1C1 Polyclonal Antibody For Wb, Elisa - Q04828

Item Number: CM0017656



Introduction

AKR1C1 rabbit polyclonal antibody validated for WB and ELISA, reacting with human and mouse. Ideal for steroid metabolism and neurosteroid research; sourced from recombinant immunogen covering full-length human AKR1C1.

[Learn More](#)

Parameter	Value
Product Name	AKR1C1 Rabbit pAb
Remarks/Alias	C9; DD1; DDH; DDH1; H-37; HBAB; MBAB; HAKRC; DD1/DD2; 2-ALPHA-HSD; 20-ALPHA-HSD; AKR1C1
Species	Human
GeneID (Human)	1645
GeneID	1645
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human AKR1C1 (NP_001344.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q04828
Protein Weight	37kDa
Shipping	Ice bag