


Anti-Akr1C2 Polyclonal Antibody (Rabbit Pab) For Wb, Ihc-P, If/Icc, Elisa - P52895

Item Number: CM0014354



Introduction

CamelBio offers high-quality anti-AKR1C2 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat; detects 37kDa protein. Essential tool for steroid metabolism and hormone research.

[Learn More](#)

Parameter	Value
Product Name	AKR1C2 Rabbit pAb
Remarks/Alias	DD; DD2; TDD; BABP; DD-2; DDH2; HBAB; HAKRD; MCDR2; SRXY8; DD/BABP; AKR1C-pseudo; AKR1C2
Species	Human
GeneID (Human)	1646
GeneID	1646
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human AKR1C2 (NP_995317.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52895
Protein Weight	37kDa
Shipping	Ice bag