

Anti-Akr1C3 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P42330

Item Number: CM0020497



Introduction

High-quality rabbit polyclonal antibody against human AKR1C3, validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with mouse. Suitable for researching androgen metabolism, steroid hormone signaling, and hormone-dependent cancers.

[Learn More](#)

Parameter	Value
Product Name	AKR1C3 Rabbit pAb
Remarks/Alias	DD3, DDX, PGFS, HAKRB, HAKRe, HA1753, HSD17B5, hluPGFS, AKR1C3
Species	Human
GeneID (Human)	8644
GeneID	8644
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42330
Protein Weight	37kDa
Shipping	Ice bag