

Anti-Akr1C3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P42330

Item Number: CM0010163



Introduction

High-quality rabbit monoclonal antibody targeting human AKR1C3 (P42330). Validated for WB, IF/ICC, ELISA; cross-reacts with human and mouse. Essential for androgen metabolism and hormone-dependent cancer research. Detects ~37kDa band.

[Learn More](#)

Parameter	Value
Product Name	AKR1C3 Rabbit mAb
Remarks/Alias	DD3; DDX; PGFS; HAKRB; HAKRe; HA1753; HSD17B5; hluPGFS; AKR1C3
Species	Human
GeneID (Human)	8644
GeneID	8644
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human AKR1C3 (P42330)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42330
Protein Weight	37kDa
Shipping	Ice bag