

Anti-Akr1C4 Rabbit Polyclonal Antibody For Wb, Elisa - P17516

Item Number: CM0030817



Introduction

AKR1C4 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, and rat AKR1C4 (P17516). Ideal for studying steroid hormone metabolism and neurosteroid pathways.

[Learn More](#)

Parameter	Value
Product Name	AKR1C4 Rabbit pAb
Remarks/Alias	C11; CDR; DD4; CHDR; DD-4; HAKRA; 3-alpha-HSD; AKR1C4
Species	Human
Gene ID (Human)	1109
Gene ID	1109
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17516
Protein Weight	37kDa
Shipping	Ice bag