


Anti-Akr7A3 Rabbit Polyclonal Antibody For Wb, Elisa - 095154

Item Number: CM0031050



Introduction

Rabbit polyclonal antibody against human AKR7A3, suitable for western blot and ELISA. Cross-reacts with mouse and rat. Targets aflatoxin B1 aldehyde reductase involved in detoxification. Ideal for studying cellular detoxification and toxin metabolism.

[Learn More](#)

Parameter	Value
Product Name	AKR7A3 Rabbit pAb
Remarks/Alias	AFAR2; AKR7A3
Species	Human
GeneID (Human)	22977
GeneID	22977
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-331 of human AKR7A3 (NP_036199.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95154
Protein Weight	37kDa
Shipping	Ice bag