

Anti-Akt1 (Phospho T308) Rabbit Mab For Wb, Elisa - P31749

Item Number: CM0001387



Introduction

Phospho-specific rabbit monoclonal antibody against human AKT1 phosphorylated at Thr308. Validated for Western blot and ELISA; ideal for PI3K/AKT signaling, apoptosis, and cell metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Akt-T308 Rabbit mAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-Akt-T308
Species	Human
GeneID	207
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T308 of human Phospho-Akt-T308 (NP_005154.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag