

Anti-Akt1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P31749

Item Number: CM0030275



Introduction

Rabbit polyclonal antibody targeting human AKT1 (PKB/RAC-alpha). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. Ideal for AKT signaling research in cell metabolism, proliferation, and survival.

[Learn More](#)

Parameter	Value
Product Name	AKT1 Rabbit pAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; AKT1
Species	Human
GeneID (Human)	207
GeneID	207
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human AKT1 (NP_005154.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag