

Anti-Akt1/2 Rabbit Monoclonal Antibody For Wb, Ihc, Elisa - P31749 / P31751

Item Number: CM0010853



Introduction

Rabbit monoclonal antibody targeting human AKT1 and AKT2, validated for WB, IHC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying AKT signaling in cancer and metabolic research.

[Learn More](#)

Parameter	Details
Product Name	AKT1/2 Rabbit mAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; PKBB; PRKBB; HIHGHH; PKBBETA; RAC-BETA; AKT1/2
Species	Human
GeneID (Human)	207/208
GeneID	207/208
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749/P31751
Protein Weight	56kDa
Shipping	Ice bag