

Anti-Akt1/Akt2/Akt3 Rabbit Polyclonal Antibody For Wb, Elisa - P31749/P31751/Q9Y243

Item Number: CM0024318



Introduction

Rabbit polyclonal antibody against AKT1/AKT2/AKT3, validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying AKT signaling in cancer, metabolism, and cell survival.

[Learn More](#)

Parameter	Value
Product Name	AKT1/AKT2/AKT3 Rabbit pAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA
Species	Human
GeneID (Human)	207/208/10000
GeneID	207/208/10000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-123 of human AKT1/AKT2/AKT3 (NP_005154.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749/P31751/Q9Y243
Protein Weight	56kDa
Shipping	Ice bag