


## Anti-Akt1S1/Pras40 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - Q96B36

Item Number: CM0030710



### Introduction

Rabbit polyclonal antibody targeting AKT1S1/PRAS40, a key mTORC1 negative regulator. Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Molecular weight 27kDa. Useful for cell growth and cancer research.

[Learn More](#)

| Parameter        | Value                       |
|------------------|-----------------------------|
| Product Name     | AKT1S1/PRAS40 Rabbit pAb    |
| Remarks/Alias    | Lobe; PRAS40; AKT1S1/PRAS40 |
| Species          | Human                       |
| GeneID (Human)   | 84335                       |
| GeneID           | 84335                       |
| Immunogen        | Synthetic peptide           |
| Source           | Rabbit                      |
| Category         | Polyclonal Antibodies       |
| Application      | WB, IHC-P, IF/ICC, ELISA    |
| Cross Reactivity | Human, Mouse, Rat           |
| SWISS            | Q96B36                      |
| Protein Weight   | 27kDa                       |
| Shipping         | Ice bag                     |