


Anti-Akt1S1/Pras40 Rabbit Polyclonal Antibody For Wb, Elisa - Q96B36

Item Number: CM0020567



Introduction

Rabbit polyclonal antibody targeting AKT1S1/PRAS40, a negative mTORC1 regulator, for Western blot and ELISA. Recognizes human AKT1S1/PRAS40; predicted molecular weight 27 kDa. Ideal for nutrient-sensing and cancer research.

[Learn More](#)

Parameter	Value
Product Name	AKT1S1/PRAS40 Rabbit pAb
Remarks/Alias	Lobe; PRAS40; AKT1S1/PRAS40
Species	Human
Gene ID (Human)	84335
Gene ID	84335
Immunogen	Synthetic peptide, a sequence within amino acids 150-256 of human AKT1S1/PRAS40 (NP_001092103.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96B36
Protein Weight	27 kDa
Shipping	Ice bag