


Anti-Akt2 Monoclonal Antibody For Wb - P31751

Item Number: CM0006921



Introduction

Rabbit monoclonal antibody targeting human AKT2 (P31751), validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for studying PI3K/AKT signaling, glucose metabolism, and cell survival in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	AKT2 Rabbit mAb
Remarks/Alias	PKBB; PRKBB; HIHGHH; PKBBETA; RAC-BETA; AKT2
Species	Human
Gene ID (Human)	208
Gene ID	208
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 390-490 of human AKT2.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31751
Protein Weight	56kDa
Shipping	Ice bag