

Anti-Akt3 Monoclonal Antibody For Wb, Elisa - Q9Y243

Item Number: CM0007622



Introduction

Rabbit monoclonal antibody targeting human AKT3, validated for WB and ELISA. Reacts with human, mouse, and rat AKT3. Ideal for studies of AKT signaling, cancer, metabolism, and brain development.

[Learn More](#)

Parameter	Value
Product Name	AKT3 Rabbit mAb
Remarks/Alias	MPPH; PKBG; MPPH2; PRKBG; STK-2; PKB-GAMMA; RAC-gamma; RAC-PK-gamma; AKT3
Species	Human
Gene ID (Human)	10000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 93-154 of human AKT3 (NP_005456.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y243
Protein Weight	56kDa
Shipping	Ice bag