

Anti-Akt3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y243

Item Number: CM0017904



Introduction

Rabbit polyclonal antibody targeting human AKT3 (RAC-gamma serine/threonine-protein kinase), validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for AKT signaling research in brain development and cancer.

[Learn More](#)

Parameter	Specification
Product Name	AKT3 Rabbit pAb
Remarks/Alias	MPPH, PKBG, MPPH2, PRKBG, STK-2, PKB-GAMMA, RAC-gamma, RAC-PK-gamma, AKT3
Species	Human
Gene ID (Human)	10000
Gene ID	10000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-190 of human AKT3 (NP_859029.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y243
Protein Weight	56kDa
Shipping	Ice bag