


Anti-Alad Rabbit Polyclonal Antibody For Wb, Elisa - P13716

Item Number: CM0030711



Introduction

High-quality rabbit polyclonal antibody targeting human ALAD (delta-aminolevulinic acid dehydratase/porphobilinogen synthase). Validated for Western blot and ELISA. Suitable for heme biosynthesis and porphyrin metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ALAD Rabbit pAb
Remarks/Alias	PBG5; ALADH; ALAD
Species	Human
GeneID (human)	210
GeneID	210
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human ALAD (NP_000022.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13716
Protein Weight	36kDa
Shipping	Ice bag