


Anti-Alad Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13716

Item Number: CM0018120



Introduction

Rabbit polyclonal antibody against human delta-aminolevulinic acid dehydratase (ALAD), validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Expected band 36 kDa. For heme biosynthesis research.

[Learn More](#)

Parameter	Specification
Product Name	ALAD Rabbit pAb
Remarks/Alias	PBGS; ALADH; ALAD
Species	Human
GeneID (human)	210
GeneID	210
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human ALAD (NP_000022.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13716
Protein Weight	36kDa
Shipping	Ice bag