

Anti-Alas1 Monoclonal Antibody For Wb, Elisa - P13196

Item Number: CM0008224

Introduction

Anti-ALAS1 rabbit monoclonal antibody validated for WB and ELISA. Detects human ALAS1, a mitochondrial enzyme catalyzing the rate-limiting step of heme biosynthesis. Ideal for heme metabolism and porphyria research.

[Learn More](#)

Parameter	Value
Product Name	ALAS1 Rabbit mAb
Remarks/Alias	ALAS; MIG4; ALAS3; ALASH; ALAS-H; ALAS1
Species	Human
GeneID (Human)	211
GeneID	211
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 541-640 of human ALAS1 (P13196).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13196
Protein Weight	71kDa
Shipping	Ice bag