

Anti-Alas1 Rabbit Polyclonal Antibody For If/Icc And Elisa - P13196

Item Number: CM0014145

Introduction

High-quality rabbit polyclonal antibody against human ALAS1, validated for IF/ICC and ELISA. Specifically detects human ALAS1, a key enzyme in heme biosynthesis. Ideal for mitochondrial studies; protein weight 71kDa.

[Learn More](#)

Parameter	Value
Product Name	ALAS1 Rabbit pAb
Remarks/Alias	ALAS; MIG4; ALAS3; ALASH; ALAS-H; ALAS1
Species	Human
GeneID (Human)	211
GeneID	211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ALAS1 (NP_954635.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P13196
Protein Weight	71kDa
Shipping	Ice bag