

Anti-Alas1 Rabbit Polyclonal Antibody For Wb, Elisa - P13196

Item Number: CM0024853

Introduction

Rabbit polyclonal antibody targeting human ALAS1 (5-aminolevulinate synthase 1). Validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for heme biosynthesis research and assay development.

[Learn More](#)

Parameter	Value
Product Name	ALAS1 Rabbit pAb
Remarks/Alias	ALAS; MIG4; ALAS3; ALASH; ALAS-H; ALAS1
Species	Human
Gene ID (Human)	211
Gene ID	211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ALAS1 (NP_954635.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13196
Protein Weight	71kDa
Shipping	Ice bag