

Anti-Alas2 Polyclonal Antibody For Wb, Elisa - P22557

Item Number: CM0026184



Introduction

Rabbit polyclonal antibody against human ALAS2, validated for WB and ELISA. Recognizes human, mouse, rat ALAS2; suitable for heme biosynthesis and erythropoiesis research.

[Learn More](#)

Parameter	Value
Product Name	ALAS2 Rabbit pAb
Remarks/Alias	ASB; ANH1; XLSA; ALASE; XLDPP; XLEPP; ALAS-E; SIDBA1; ALAS2
Species	Human
GeneID (Human)	212
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-320 of human ALAS2 (NP_001033057.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22557
Protein Weight	65kDa
Shipping	Ice bag