


Anti-Albumin Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P02768

Item Number: CM0018651



Introduction

Albumin Rabbit Polyclonal Antibody validated for WB, IF/ICC, IF-P, ELISA. Reacts with human, mouse, rat. Immunogen: recombinant protein. Ideal for plasma osmotic pressure and metal transport studies.

[Learn More](#)

Parameter	Value
Product Name	Albumin Rabbit pAb
Remarks/Alias	HSA; FDAHT; PRO0883; PRO0903; PRO1341; Albumin
Species	Human
GeneID (Human)	213
GeneID	213
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-609 of human Albumin (NP_000468.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02768
Protein Weight	69kDa
Shipping	Ice bag