


Anti-Albumin/Alb Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02768

Item Number: CM0005521



Introduction

Rabbit monoclonal antibody against human Albumin/ALB, validated for Western blot, immunohistochemistry, immunofluorescence, and ELISA. Suitable for human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	Albumin/ALB Rabbit mAb
Remarks/Alias	ALB; HSA; PRO0883; PRO0903; PRO1341; albumin; Albumin/ALB
Species	Human
Gene ID (Human)	213
Gene ID	213
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-609 of human Albumin/ALB (NP_000468.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02768
Protein Weight	45kDa/47kDa/69kDa
Shipping	Ice bag