

Anti-Alcohol Dehydrogenase 1A Monoclonal Antibody For Wb, Elisa - P07327

Item Number: CM0003381



Introduction

Rabbit monoclonal antibody against human Alcohol Dehydrogenase 1A (ADH1A), validated for WB and ELISA. Cross-reacts with mouse. Ideal for alcohol metabolism research. - P07327.

[Learn More](#)

Parameter	Value
Product Name	Alcohol Dehydrogenase Rabbit mAb
Remarks/Alias	ADH1; Alcohol Dehydrogenase
Species	Human
Gene ID (Human)	124
Gene ID	124
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Alcohol Dehydrogenase (P07327).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P07327
Protein Weight	40 kDa
Shipping	Ice bag