

Anti-Alcohol Dehydrogenase Rabbit Polyclonal Antibody For Wb, Elisa - P07327

Item Number: CM0014450



Introduction

Rabbit polyclonal antibody against human Alcohol Dehydrogenase 1A (ADH1A), validated for WB and ELISA. Cross-reacts with rat. Ideal for studying alcohol metabolism and ethanol oxidation.

[Learn More](#)

Parameter	Value
Product Name	Alcohol Dehydrogenase Rabbit pAb
Remarks/Alias	ADH1A; ADH1; Alcohol Dehydrogenase
Species	Human
GeneID (Human)	124
GeneID	124
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 220-300 of human Alcohol Dehydrogenase (NP_000658.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P07327
Protein Weight	40kDa
Shipping	Ice bag